



Description

Vomit support
Designed for Hygienic Covers®
Bedpan washer compatible
Recyclable
Made in Canada

Item Number

PG1016

Presentation

Polypropylene vomit support.
 It is used in association with a Hygienic cover® – vomit bag containing a super absorbent pad. It allows vomit collection.

Colour

White.

Support dimensions

Length: 7.5" (19 cm) x Width: 4.5" (11.5 cm) x Height: 5.1" (13 cm).

Composition

Polypropylene. No latex and no phthalates.

Shelf life

No expiry date.

Manufacturing process

Injection moulding

Use

For any person, men or women, experiencing vomiting, whether in a health care institution, at home or during a trip. Particularly recommended for oncology services, dialysis, post-operative recovery room, as well as for emergency services, such as ambulances, etc.

Instructions for use

Refer to the tutorial for more detail.

Precautions for use

Single-patient use. It is recommended to dispose of vomit support after the patient's departure.

Cleaning / Sterilization

The vomit support must be regularly cleaned and when soiled. It is recommended to use a disinfectant detergent or an enzyme product associated with a mechanical rubbing action. Follow the instructions on disinfectants to maximize their efficiency (dilution, contact time ...). It can also be cleaned in a bedpan washer without exceeding the temperature of 203°F or 95°C.

Disposal

Dispose of product in regular waste or according to the protocol of the institution. The support is recyclable (Plastic #5 - PP).

Benefits

Increases hygiene, saves time, reduces odors, simple and risk-free disposal.

Identification

Each case of vomit supports is identified with a lot number.

Regulatory

FDA: class I medical device exempt from 501(k).
CE marking: Class I medical device according to the rules of annex IX of European Council Directive 93/42/EEC.

Risks / Safety

This product was evaluated and is not considered hazardous based on the hazard criteria established by the OSHA (Occupational Safety and Health Administration). Safety data sheet is available upon request.

Handling/Storage

Store in a clean and dry place, away from sun.

Packaging

Unit	Box of 20 Hygienic Covers	19 X 13 X 11,5 cm / 75.4 g 7.5 x 4.5 x 5.1 in / 0.17 lbs
Case	20 vomit supports per case	27 x 60 x 22 cm / 2.3 kg 10,5 x 23,5 x 8,5 in / 5 lbs
Pallet	40 cases or 800 vomit supports	122 x 107 x 120 cm / 113 kg 48 x 42 x 47 ¼ in / 250 lbs

