

**Description** Graduated vomit bag with super-absorbent pad  
**Absorption capacity: 300 ml**  
**No expiration date**  
**Eco-friendly**  
**Assembled in Canada**



**Item Number** PG1022

**Utility** Single use. Transparent vomit bag graduated in ml with cardboard cone with non-return valve and containing a super-absorbent pad. Allows to solidify biological liquids and to reduce odors. Collects vomit.

**Envelop dimensions** Length: 39.5 cm (15.5 in.). Width: 15 cm (6 in.). Thickness: 38µm

**Composition** Clear polyethylene bag with cardboard cone. 100% recycled super-absorbent pad composed of a non-woven sulfate-based backing coupled with sodium acrylate copolymers. Packaging made of recycled materials. Does not contain latex or phthalates.

**Product Lifespan** No expiration date.

**Efficiency** Absorption of 300ml of liquid in 30 seconds.

**Manufacturing process** Thermo welding of plastic film.

**Usage** For any person experiencing vomiting either while on-the-go, traveling, in an institution, public area or at home. Particularly recommended for oncology, dialysis, post-operative recovery room or for all emergency services such as ambulances.

- Instructions for use**
1. Open the box and take out a bag.
  2. Open the cone.
  3. Vomit into the bag.
  4. Wait for the liquids to solidify (±300 ml in a few seconds).
  5. Dispose of in the appropriate place.
- Refer to the instruction manual for more details.

**Precautions for use** Single use product. Do not reuse after use. Do not leave within the reach of children.

**Disposal** Dispose of in regular garbage or according to the facility protocol. Do not flush down the toilet.

**Advantages** Improves hygiene, saves time, reduces odor, is easily and safely disposed of. Replaces the bean and paper bag. Limits the spread of pathogens. Single use product.

**Identification** Each product is identified by a sku number.

**Quality** Our bags are tested according to the methodology and acceptance criteria of NF-X30-501 standard: Packaging for healthcare waste- Bags for soft healthcare waste – Specifications and test methods.

**Regulations** **FDA:** Class I medical device exempt from 501 (k)  
**CE Marketing:** This product belongs to Class I according to the rules of Annex IX of the European Medical Device Directive 93/94/EEC.

**Risk / Security** Very slight risk of eye or respiratory tract irritation in case of prolonged contact with the super absorbent pad. Be careful not to use the pad as a wipe on the skin or in diapers (baby or incontinence). The Material Safety Date Sheet is available upon request.

**Handling / Storage** Keep clean, away from the sunlight and moisture.

**Packaging**

Unit	Box of 50 vomit bags	19 x 20 x 14,5 cm / 1,3 kg 7,5 x 8 x 5,75 po / 2,9 lb
Case	8 boxes of 50 vomit bags per case.	60 x 27 x 42 cm / 12,3 kg 23,5 x 10,75 x 16,5 po / 27,2 lb
Pallet	20 cases or 160 boxes of 50 vomit bags.	126 x 102 x 121 cm / 269 kg 49,5 x 40 x 47,75 pouces / 594 lb

